

Care Under Pressure

ReddyPort™ non-invasive ventilation products

Non-invasive ventilation (NIV) is the first line of therapy in respiratory insufficiency or failure, commonly seen with COPD, CHF, Asthma, pneumonia or ARDS¹, yet many patients struggle because the mask causes dry-mouth, phlegm build up, and an inability to communicate.

Patient discomfort can cause NIV failure, which leads to mechanical ventilation and the escalation of care, increasing risks and costs.

ReddyPort products promote NIV success, by empowering clinicians and patients with simplified oral access while maintaining pressure. This access allows new care to be provided while therapeutic pressure is maintained, improving patient compliance and comfort, leading to more successful NIV.



ReddyPort 531 elbow



ReddyPort 541 elbow

ReddyPort™ Elbow

The ReddyPort Elbow replaces the standard elbow on an NIV mask. Featuring a patented self-sealing valve, the ReddyPort elbow fits onto compatible masks, providing easy access to the oral cavity while maintaining therapeutic pressures.

- Patented self-sealing valve seals with ventilator pressure and self-reverts
- Includes an anti-asphyxia valve and connects to standard NIV tubing
- ReddyPort RP531 elbow connects to Philips Respironics AF531 and PerfoMax NIV masks
- ReddyPort 541 elbow connects to the Philips Respironics AF541 NIV mask

"Being able to clean the mouth is such a blessing. She cannot be removed from the mask because of shortness of breath. We are very happy."

Family of patient at Huntsman Cancer Hospital. Adult female with respiratory failure from pneumonia pulmonary edema, DNI status.



ReddyPort toothbrush being used provide oral hygiene care



The kit includes swabs, toothbrushes, a yankauer (featured in photo above) oral rinses, and moisturizing gel (not shown).

ReddyPort™ NIV Maintenance Kit

The ReddyPort NIV Maintenance Kit includes everything needed for oral hygiene and keeping the patient's mouth clean and moisturized, without the risk and time required to remove the patient's mask.

- Patented components seal with ReddyPort Elbow
- NIV pressure is maintained while cleaning and moisturizing, reducing mask removal
- Protective guard allows for easy cleaning and moisturizing, enabling appliance insertion and removal
- Decreases dependency on clinicians for issues like dry mouth
- Addresses pneumonia risk factors by improving compliance to oral care protocols
- The NIV Maintenance Kit follows Q4 protocols to allow for standard oral care procedures to be applied to NIV patients



Patient using microphone to communicate



Microphone with plug-in speaker and Controller

For more information about ReddyPort products
call **801.899.3036**.

Reference

1 Rochwerg B, Brochard L, Elliott MW, Hess D, et al. Official ERS/ATS Clinical Practice Guidelines: Noninvasive Ventilation for Acute Respiratory Failure. Eur Respir J. 2017; 31;50(2):1602426.

ReddyPort, Salt Lake City, UT, 84101, U.S.

reddyport.com

© 2020 ReddyPort and the ReddyPort logo are trademarks of ReddyPort and Company. FDA cleared, patented self-sealing elbow, oral care and other patents pending. The Philips and Respironics trademarks are property of Koninklijke Philips N.V. or their respective owners. RP2000 (0420)

REDDYPORT